



huka

TAKE A
FREE
TESTRIDE!

TRICYCLE CITY

The sporty Dutch design city bike

A lightweight bicycle with stability

The City tricycle features a new contemporary look from Dutch Design with even more cycling comfort. The bicycle has been developed for users who need more stability and extra safety for frequent use. The low entry height ensures that you can easily get on and off. The tricycle can also be adjusted to your wishes and circumstances. There are various accessories available with which you can make your bike even more personal. The City makes you more mobile and independent!

We look at individual cyclists and their needs and want to match them as well as possible. This way, we continuously innovate and optimize our bikes, so that your bike is equipped with the best materials and optimum comfort. For example, the City is light in weight and also suitable for taller cyclists, due to the increased crank length.

The City tricycle comes standard with a parking brake, a rear drum brake, a largely closed modern chain guard and a safety lock (AXA solid).

- ✔ Available with electrical support and available in two sizes.
- ✔ Extremely stable and manoeuvrable due to the low centre gravity.
- ✔ The adjustment range of the saddle is extended, so that an optimal sitting position can be achieved.
- ✔ 10-year warranty on the frame and the front fork.*

* Register your Huka bike once via www.huka.nl/en/warranty and get a 10 year warranty on the frame and front fork. Standard you get a 5 year warranty on these parts and 2 years on other components (with the exception of wear sensitive parts).

Electrical support

The City is available with pedal support (PAS-vario). With the help of this drive it is easier to pedal and therefore less tiring. This silent pedal support is standard supplied with a start help, five travel speed settings and a clear display.

In addition, the City is also workable with AUTO-Vario. This allows you to ride without pedalling. The cycling speed is measured via the gas handle.

Support & special adjustments

With the City we can offer various extra supports such as a backrest, various foot fixations, etc. This way we can make specific adjustments to your bike.

Accessories & options

The City can be equipped with various accessories and options. Do you want a coaster brake, a reverse pedalling automatic or an 8-speed hub? View the options on our website or inquire with your dealer.

Colours

The tricycle City standard comes in satin light grey (RAL7040).



The 10 most important specifications

Technical specifications	City, size M	City, size L
Length in mm	1675	1675
Width in mm	785	785
Weight in kg	25	26
Crank length in mm	170	170
Length inner leg in mm	660 - 880	790 - 980
Entry frame length in mm	570	570
Entry height frame in mm	190	190
Maximum load in kg	130	130
Tire size front in inch / ETRTO	24 x 1,75/ 47-507	24 x 1,75/ 47-507
Tire size rear in inch / ETRTO	20 x 1,75/ 47-406	20 x 1,75/ 47-406

About Huka

Developed in Oldenzaal, we have been producing bikes that make a difference in the world of riders for more than 40 years. Your mobility, independence and freedom is our passion. With a committed team, we provide a wide collection of innovative bicycles and other means of transport with Dutch design. You want to learn more about the City and our assortment? Give us a call or visit our website!

