



**huka**

MAKE A  
FREE  
TESTRIDE!

**WHEELCHAIR SCOOTER PENDEL**

# The scooter for wheelchair users

**Independently, safely & comfortably on the road**

Are you dependent on a wheelchair and would you like to cover greater distances independently and safely? Meet the Pendel wheelchair scooter! The wheelchair scooter is the alternative to the car or (regional) taxi. The Pendel is easy to use, because you can drive around completely independently with your own wheelchair! A strenuous switch is therefore not necessary! Do you have additional wishes tailored to your circumstances? Various options are possible! With a maximum speed of 25 km/h and a radius of action of up to 60 kilometres, you will reach your destination quickly and safely. Such freedom!

We look at individual drivers and their needs and want to match them as well as possible. This way, we continuously innovate and optimize our wheelchair scooters, so that your wheelchair scooter is equipped with the best materials and optimum comfort.

The Pendel is equipped with a large battery capacity, powerful mechanical and hydraulic brakes and a dashboard with display. The large headlight offers an enormous light output, so that you have a good view while driving and you are also well seen by other road users.

- ✔ Easy and fast travel with the wheelchair.
- ✔ The suspension fork and the air suspension at the rear ensure comfortable driving.
- ✔ The electric drive makes the Pendel quiet, emission-free and inexpensive to use.
- ✔ Complies with strict European safety requirements and has been approved as a vehicle and medical device.

## Electrical support

The Pendel is standard equipped with electrical support. Because of the mechanical and hydraulic braking system you are able to brake in a controlled and powerful way. This provides a short braking distance.

## Support & special adjustments

Do you have specific needs to be able to ride optimally? We can make specific individual adjustments to your wheelchair scooter.

Your special need is a challenge for us!

## Accessories & options

The Pendel can be styled entirely to your liking and taste with various accessories and options. Do you want an extra mirror, transfer handles or a custom-made ramp entirely tailored to your wheelchair? Various options are possible. That is no problem for us! View the options on our website or inquire with your dealer.

## Colours

The Pendel wheelchair scooter comes in 6 standard colours:

- Black (RAL 9017)
- Blue (RAL5025)
- Grey (RAL7043)
- Green (RAL6028)
- Red (RAL3001)
- Orange (Y4-Y436)

Would you like a different RAL colour? This is available on request.



## The 10 most important specifications

Technical specifications	Pendel, standard version
Maximum wheelchair width in mm	800
Total length in mm	1760
Total width in mm	1000
Turning circle in mm	3900
Maximum rise in degrees	170
Maximum load in kg	125
Total weight Pendel incl. batteries in kg	160
Tire size front/ rear	3.00-10/ 3.00-12
Maximum speed as a medical device/ vehicle in km/h	15/ 25
Range in km	ca. 40-60

## About Huka

Developed in Oldenzaal, we have been producing bikes that make a difference in the world of riders for more than 40 years. Your mobility, independence and freedom is our passion. With a committed team, we provide a wide collection of innovative bicycles and other means of transport with Dutch design. You want to learn more about the Pendel and our assortment? Give us a call or visit our website!

