





Technical data AT-B

	Standard version
Material frame	anti-corrosion treated steel
Coating	outside quality epoxy powder coating
Color ^A	blue (RAL5015) or pink (RAL3014)
Wheels	aluminium rim and galvanized hub with stainless steel spokes
Tires	Schwalbe
Lighting	LED battery lighting front and rear
Bell ^A	sporty and easy-to-use clock
Mudguards	impact and form-retaining, plastic covered
Chain protection	impact-resistant closed plastic chain guard
Handlebar and stem ^A	slightly bent, adjustable, stem 300 mm
Handles	rubber
Seat ^A	children's seat, Selle Royal, height adjustable
Seatpost ^A	aluminium, 350 mm, height adjustable
Brakes	roller brake: front- and rear, parking brake
Drive ^A	freewheel with differential in rear axle
Key lock	AXA Solid safety lock
Options	see price list for different options

	Standard version
Total length in mm (L)	1360
Width total in mm	625
Weight in kg	24
Crank length in mm	110
Length inner leg in mm (X)	470-650
Seat height to ground in mm (Z)	min. 530
Entry space in mm (Y)	340
Entry height in mm (H)	260
Height adjustment handlebar in mm	70
Handlebar width in mm	570
Maximum load in kg	60
Tire size front ETRTO	47-406
Tire size front in inches	20 x 1,75
Tire size rear in ETRTO	47-305
Tire size rear in inches	16 x 1,75



Warranty

Huka has a standard 5-year warranty on the frame and front fork. If you register your bike via www.huka.nl/en/warranty, you are eligible to receive a 10-year warranty on these parts. Additionally you have a 2-year warranty on other parts with the exception of wear-sensitive parts. You can find more information about the AT-B warranty on www.huka.nl/en/warranty.

Possibilities for alternative parts

The AT-B is equipped with a number of standard parts. Optional these parts can be replaced by alternative parts. The relevant parts are marked by a symbol behind this specific part. The alternative parts can be found on the price list

The delivered product may differ from the images shown. The accessories shown are not included. Printing errors reserved.