

Pricelist/ order form City



Name _____

Organisation _____

Street/ Nr. _____

Postal code/ City _____

Telephone _____

Base City

The City basically comes with a number of standard parts. These parts are marked with a symbol in front of this specific part. The City will be provided with these parts if no other option from this category has been selected. Also a number of parts come with a standard addition. These parts are marked with a symbol behind. It will have this standard addition if no alternative has been selected. For the complete description of the standard parts we refer to the specifications. Below you will also find all possible options and accessories.

Art.nr.	Base	Price €
401600	Tricycle City, size M	1.989
402620	Ⓢ - Frame M, color grey (RAL7040)	-
402630	- Custom RAL color frame M (excl. fluorescent colors), additional charge	199
401700	Tricycle City, size L	2.132
402720	Ⓢ - Frame L, color grey (RAL7040)	-
402730	- Custom RAL color frame L (excl. fluorescent colors), additional charge	199

Art.nr.	Drive options/ hubs	Price €
Ⓢ 409400	Freewheel	-
409401	Fixed gear	-
409402	Pedal brake	129
409403	Pedal brake automatic	283
409404	8 speed hub freewheel (Shimano Nexus)	222
409405	8 speed hub with pedal brake (Shimano Nexus)	299
409406	8 speed hub with pedal brake automatic (Shimano Nexus)	427
409407	Ⓢ - Gripshift for serving 8-speed hub (Shimano Nexus)	-
409408	Thumb adjuster for serving 8 speed hub (Shimano Nexus)	83

Art.nr.	(Electrical) support	Price €
Ⓢ 409420	Without electrical support	-
401901W	Pedal support PAS-vario Bafang 43 Volt (pedaling is necessary)	2.515
401930	Ⓢ - Base accupack 450Wh (range approx. 40-50 km)	-
401932	- Upgrade battery pack to 750Wh (range approx. 60-70 km)	255

Art.nr.	Braking	Price €
Ⓢ 409410	Rollerbrake at rear wheel (applicable with freewheel of 8 speed hub freewheel (Shimano)*)	-
412838L	- Brake lever for double braking, left**	93
412838R	- Brake lever for double braking, right**	93
406253	Hydraulic rear disc brakes, control right side	252
406251	Hydraulic rear disc brakes, control left side	252

* This braking is not possible in combination with rear hydraulic disc brakes.

** This braking is not possible in combination with a pedal brake (automatic) and/ or rear hydraulic disc brakes.

Make the tricycle City even more personal? Special adjustments are possible on request!

The mentioned prices are exclusive VAT, delivery ex works. The prices are valid from January 1, 2020 and replace all previous price lists. Subject to technical changes and printing errors.

Pricelist/ order form City

	Art.nr.	Handlebars	Price €
S	401424	Sporty slightly curved handlebar, black	-
	409205	Handlebar butterfly model	61
	406148	Handlebar, rectangular model	145
	430726	Handlebar, rectangular model incl. push braking system*	381

* not possible in combination with rear hydraulic disc brakes

	Art.nr.	Stems	Price €
S	409039	Ergonomic stem	-
	407120	Ergonomic stem with quick release	92

	Art.nr.	Other handlebar attachments	Price €
S	426827	Bicycle clock type Decibel*	-
	426959	Rotation clock	12
	427400	Handlebar pointed backwards (steer closer)	58
	426958	Extra stroke limitation, for handlebar	22
	426961	Rearview mirror	32
	426924	Lower arm support	161
	426927	Lower arm fixation (only suitable for the item above)	33
	426951	Support bracket, handlebar	100
	426992	Vertical handle, black	39

*This clock is always delivered as standard, outside the other selected parts in this category. Unless a rotation clock is chosen.

	Art.nr.	Seats	Price €
S	401637	Velo Plush seat	-
	426833	Female seat, synthetic version	-
	426837	Female seat, leather version	147
	426836	Male seat, leather version	147
	427179	Gel seat	58
	418265	Wide seat	58
	426897	Banana seat (not suitable with T-model seat post)	41
	418266	Puch seat	84

	Art.nr.	Seat posts	Price €
S	409033	Aluminium seat post, length 350 mm	-
	432059	Suspension seat post, length 350 mm	48
	401591	T-model seat post, length 400 mm	63
	401596	Extended and reinforced seat post, length 400 mm	67

	Art.nr.	Trunk support	Price €
	400011	Backrest with hip belt (adjustable in height and width)*	210
	409012	S - Seatpost for the purpose of adapter back support**	54
	400013	Back support and lumbar support with belly belt (adjustable in height and width) *	267
	409012	S - Seatpost for the purpose of adapter back support**	54

* Pay attention! The required inner leg length will be 5 centimeters longer when this support is used. With this support the seatpost for the adapter of the backrest have to be ordered.

** The selected seat post in the seat post category is hereby cancelled.

Make the tricycle City even more personal? Special adjustments are possible on request!

The mentioned prices are exclusive VAT, delivery ex works. The prices are valid from January 1, 2020 and replace all previous price lists. Subject to technical changes and printing errors.

Pricelist/ order form City

Art.nr.	Stick holders	Price €
406155	Stick holder, double rear	98
406156	Stick holder, single rear	89

Art.nr.	Pedals	Price €
S 401645	Pedal set	-
426994	Balance pedals (not suitable with footrest/ forefoot support), per pair	35
426993	Large pedals (7.5 cm deep, 10 cm wide)	18
426929	Forefoot support, left (toeclip)	22
426930	Forefoot support, right (toeclip)	22

Art.nr.	Foot support	Price €
426938	Plastic footrest, left with forefoot- and heel strap (22 x 8 cm) (applicable with basic pedal set)	103
426939	Plastic footrest, right with forefoot- and heel strap (22 x 8 cm) (applicable with basic pedal set)	103
426940	Plastic footrest, left with forefoot- and heel strap (24 x 9 cm) (applicable with basic pedal set)	103
426941	Plastic footrest, right with forefoot- and heel strap (24 x 9 cm) (applicable with basic pedal set)	103
426986	Pedal extender, left	53
426987	Pedal extender, right	53
432750	Balance weight for footrest/ forefoot support, left (not suitable with lower leg fixation)	83
432751	Balance weight for footrest/ forefoot support, right (not suitable with lower leg fixation)	83
426942	Adjustable footrest metal, unbreakable on the left (applicable with pedal shaft extender)	113
426943	Adjustable footrest metal, unbreakable on the right (applicable with pedal shaft extender)	113
406174	Crank shortener set (reduces inner leg length to 2,5 cm or 5 cm)	79
428931	Metal footrest for lower leg fixation, left (25 x 9 cm) (applicable with basic pedal set)	110
428923	Lower leg fixation, left (only suitable for the above item)	79
428930	Metal footrest for lower leg fixation, right (25 x 9 cm) (applicable with basic pedal set)	110
426949	Lower leg fixation, right (only suitable for the above item)	79
432789	Pulley system/ pointed foot correction (applicable with footrests)	71

Art.nr.	Baskets & bags	Price €
427025	Shopping basket incl. clamp on handlebar	73
496002	Double bicycle bag Tour, black	54

Art.nr.	Other adjustments	Price €
426578	Safety lock (served single-handed)	22
426587	AXA Cable lock (served single-handed)	34
496055	Direction indicator set (mounted at the rear)	97
401675	Adapter for walker holder	144

Date of order

Signature/ stamp

Total €

VAT

Total € incl. VAT

Make the tricycle City even more personal? Special adjustments are possible on request!

The mentioned prices are exclusive VAT, delivery ex works. The prices are valid from January 1, 2020 and replace all previous price lists. Subject to technical changes and printing errors.