





## Technical data Duet

	Standard version
Material frame	anti-corrosion treated steel
Coating	outside quality epoxy powder coating
Wheels	aluminium rim and hub with stainless steel spokes
Tires	Schwalbe Marathon Plus (anti-puncture)
Lighting	LED battery lighting front and rear
Bell	sporty and easy-to-use bell
Mudguards	impact and form-retaining, plastic covered
Chain protection	impact-resistant closed plastic chain guard
Handlebar & seat post <sup>A</sup>	metal, suspension seat post 350 mm
Seat <sup>A</sup>	sporty gel seat, height adjustable
Brakes	drum brake front, coaster brake rear
Hub <sup>A</sup>	3-speed hub with coaster brake
Bike part	dual leg stand, seat post and cable lock
Wheelchair part*	seat shell with belt and footrest
Disconnectable frame	standard with quick coupling
Key lock	AXA cable lock
Options	see price list for different options

\* Including height adjustable and removable footboard, washable padding and pelvic band.

	Bicycle part	Wheelchair part
Total length in mm (H)	2350	
Width total in mm (M)	650	
Weight in kg	41	
Maximum load total in kg	180	
Maximum load in kg	100	100
Crank length in mm	170	
Length inner leg in mm (F) (D)	780-930	100-550
Seat height to ground in mm	910-1060	
Entry depth in mm	520	
Entry height in mm	390	
Handlebar height in mm (G)	810-1110	
Seat width in mm (A)		450
Seat depth in mm (B)		400
Back height in mm (C)		550
Tire size in ETRTO	47-559	47-559
Tire size in inch	26 x 1,75	26 x 1,75
Size front wheels		200 x 50



### Warranty

Huka has a standard 5-year warranty on the frame and front fork. If you register your bike via [www.huka.nl/en/warranty](http://www.huka.nl/en/warranty), you are eligible to receive a 10-year warranty on these parts. Additionally you have a 2-year warranty on other parts with the exception of wear-sensitive parts. You can find more information about the Duet warranty on [www.huka.nl/en/warranty](http://www.huka.nl/en/warranty).

### Possibilities for alternative parts

The Duet is equipped with a number of standard parts. Optional these parts can be replaced by alternative parts. The relevant parts are marked by a symbol behind this specific part. The alternative parts can be found on the price list.

The delivered product may differ from the images shown. The accessories shown are not included. Printing errors reserved.