

The Huka logo is a yellow diagonal banner with the word 'huka' in a dark, lowercase, sans-serif font.A yellow circular badge with the text 'TAKE A FREE TESTRIDE!' in a bold, black, sans-serif font.

**CUSTOM CHILDREN AT-SERIES**

**Happy & safe for children**

## Specialty developed tricycle for children

Are you looking for a children's bike that offers more stability or support? With the adapted AT bike your child can enjoy cycling without difficulties. Playing outside, cycling to school and playing with friends? Easy! This children's bicycle can also be adjusted entirely to the wishes and circumstances of your child. For example, the lightweight bicycle can be supplied in all hip, cool and cheerful RAL colours and there are various accessories with which you can make the bicycle even more personal and attractive. The children can cycle with pleasure and don't have to think about anything else!

We look at individual cyclists and their needs and want to match them as well as possible. This way, we continuously innovate and optimize our bikes, so that your bike is equipped with the best materials and optimum comfort. We have two types that have been specially developed for two age groups, the AT-A (up to and including about 5 years) and the AT-B bike (from about 5 years up to and including 9 years).

The AT bike is standard equipped with, among other things, a closed modern chain guard, lightly operated Rollerbrake hand brakes and a low and wide step-on plate. The AT-B is also equipped with lighting and a safety lock (AXA solid).

- ✔ Safely getting on and off, due to the three-wheel construction and the low centre gravity point.
- ✔ Smooth driving bicycle with a stable road holding.
- ✔ Individually adaptable with various supports and accessories.
- ✔ Cheerful and child-friendly styled.
- ✔ 10 year warranty on the frame and front fork.\*

\* Register your Huka bike once via [www.huka.nl/en/warranty](http://www.huka.nl/en/warranty) and get a 10 year warranty on the frame and front fork. Standard you get a 5 year warranty on these parts and 2 years on other components (with the exception of wear sensitive parts).



### Electrical support

The ATB is available with an auxiliary motor (PASvario). With the help of this drive it is easier to pedal and therefore less tiring. This silent pedal support is supplied with three travel speed settings and a well-arranged display.

### Support & special adjustments

In addition to the extra support that we can offer, such as a push bar, hip support, trunk support, a back support, foot

fixations, etc. We also provide customized solutions. This allows us to make specific adjustments to the children's bicycle. Special needs are a challenge for us!

### Accessories & options

The AT can be equipped with various accessories and options. Think of a different seat, handlebars or pedals. View the options on our website or inquire with your supplier.

### Colours

The AT-A tricycle is standard supplied in light blue (RAL5024). The AT-B tricycle standard comes in the following colours:

- Bright blue (RAL5015)
- Pink (RAL3014)

Want to make the AT even more personal with, for example, your child's favourite colours? That is possible! We can make the AT children's bike in any hip, cool or cheerful RAL colour.



### The 10 most important specifications

Technical specifications	AT-A	AT-B
Length in mm	1140	1360
Width in mm	540	625
Weight in kg	17	24
Crank length in mm	90	110
Length inner leg in mm	380-490	470-650
Entry frame length in mm	300	340
Entry height frame in mm	145	260
Maximum load in kg	40	60
Tire size front in inch / ETRTO	16 x 1,75/ 47-305	20x 1,75/ 35-406
Tire size rear in inch / ETRTO	12 x 1,95/ 50-203	16 x 1,75/ 47-305

### About Huka

Developed in Oldenzaal, we have been producing bikes that make a difference in the world of riders for more than 40 years. Your mobility, independence and freedom is our passion. With a committed team, we provide a wide collection of innovative bicycles and other means of transport with Dutch design. You want to learn more about the AT bike and our assortment? Give us a call or visit our website!