

The Huka logo is written in a bold, lowercase, sans-serif font in a dark grey color, set against a bright yellow background that is part of a diagonal banner.A circular yellow callout with the text "TAKE A FREE TESTRIDE!" in a bold, uppercase, sans-serif font.

DUO-BIKE ORTHROS

The social bike

Comfortably & independently together on the road

Independent, stable & safe on the road on a bicycle? Do you like to go out together, but you can't shake the fear of losing control of the bike? Then choose Orthros. This duobike makes it possible to be mobile and enjoy cycling, unlimited. On the Orthros, the guide and passenger are seated right next to each other. A nice chat is like that so much easier! Depending on your wishes, the Orthros will be tailored to your personal circumstances.

We look at individual riders and their needs and want to create a better world for them. We are innovating and optimizing our bikes continuously, so that your bike is equipped with the best materials and ensures you optimal comfort. It's super easy to get on the Orthros and the duo bike ensures a stable ride because of its ingenious design. This duo bike is standard equipped with a parking brake, easy to use hand brakes, hydraulic disc brakes, 8 speed freewheel, tough Big Ben anti-puncture tires and a safety lock (AXA solid).

- ✔ Standard equipped with a rotatable co-driver's seat.
- ✔ It's possible to let the passenger pedal actively, passively or not at all.

- ✔ Eye contact & a pleasant chat is easy.
- ✔ The handlebars are adjustable in height and angle to create an optimal cycling position.

A circular yellow callout with the text "NEW! COUGH SCREEN" in a bold, uppercase, sans-serif font.

Do you want to limit the risks and cycle safer on the duobike in the social distance society? Use a cough screen to protect each other.

The accessory minimizes the risk of dripping moisture from one person to another. The cough screen is made of strong plexiglass and is placed between the rider and co-driver's seat.

Electrical support

The Orthros is available with a (PAS-vario) auxiliary engine, which makes it easier to pedal and therefore it's less tiring. This silent pedal assistance is standard delivered with a drive help, three travel speed settings and one neat display.

Support & custom made

In addition to the extra supports that we can offer, as a stable footboard, armrest set etc., we deliver customization.

Therefore we are able to make specific adjustments to your bike. Your special needs are a challenge for us!

Accessories & options

The Orthros can be expanded with various accessories and options. Do you prefer another steering wheel or pedals or do you want that the passenger doesn't have to pedal at all? No problem, we will take care of your wishes!

Check out the possibilities via our website or via your dealer

Colors

The Orthros duo-bike will be standard delivered in the following colors:

- black matt (RAL9005)
- white silky gloss (RAL9016)

You want to design your Orthros even more personal? Don't worry, we got you! We can make the Orthros in any desired RAL color for you!



The 10 most important specifications

Technical specifications	Orthros
Crank length in mm	150
Inner leg length driver in mm	650-850
Inner leg length passenger in mm	650-850
Seat height (from the ground) in mm	min. 660
Tire size front in inch/ ETRTO	24 x 2,15/ 57-507
Tyre size rear in inch/ ETRTO	20 x 2,15/ 55-406
Driver maximum load in kg	120
Co-driver maximum load in kg	120
Seat width in mm per side	450
Seat depth in mm per side	380

About Huka

Developed in Oldenzaal, we have been producing bikes that make a difference in the world of riders for more than 40 years. Your mobility, independence and freedom is our passion. With a committed team, we provide a wide collection of innovative bicycles and other means of transport with Dutch design. You want to learn more about the Orthros and our assortment? Give us a call or visit our website!