

The Huka logo is written in a bold, lowercase, sans-serif font. It is white and positioned diagonally within a yellow parallelogram shape in the top-left corner of the advertisement.A yellow circular badge with the text "TAKE A FREE TESTRIDE!" in a bold, uppercase, sans-serif font. The badge is located in the top-right corner of the image, overlapping the background scene.

TANDEM COPILOT 3

Stable getting on and off with three wheels

Extra safe and stable with the passenger in front

Are you looking for a tandem that is extra stable and safe for, for example, a restless mobile passenger or a passenger with a reduced sense of balance? Experience the tandem Copilot 3, suitable for adults and children! Thanks to the third wheel, both riders can remain seated even when the tandem is standing still. With this tandem, only the supervisor can steer, and the passenger has a handlebar for guidance. The Copilot 3 can also be personalized to specific preferences and wishes. It is a tandem to enjoy cycling together without worrying about safety!

We look at individual cyclists and their needs and want to match them as well as possible. This way, we continuously innovate and optimize our bikes, so that your bike is equipped with the best materials and optimum comfort. For example, the supervisor is positioned higher than the passenger to keep a good overview of the traffic. In addition, you can easily get on and off with the Copilot 3. The supervisor on the back can optionally decide whether the passenger has to pedal along or not.

The tandem Copilot 3 is standard equipped with a low step-up entry, easy-to-use hand brakes, a large closed modern chain guard and an AXA cable lock.

- ✔ Extra stable due to two front wheels. Also, for restless moving passengers.
- ✔ Available in two models: supervisor adults & supervisor child.
- ✔ Standard equipped with a removable and tiltable front frame, making it easier to get on and to store the bike more compactly.
- ✔ 10-year warranty on the frame and the front fork.*

* Register your Huka bike once via www.huka.nl/en/warranty and get a 10 year warranty on the frame and front fork. Standard you get a 5 year warranty on these parts and 2 years on other components (with the exception of wear sensitive parts).

Electrical support

The Copilot 3 is available with an auxiliary motor (PASvario). With the help of this drive it is easier to pedal and therefore less tiring. This silent pedal support is standard supplied with a start help, three travel speed settings and a clear display.

Support & special adjustments

In addition to the extra accessories that we can offer, such as hip support, trunk support, foot fixations, etc. We also provide customized

solutions. That way we can make specific adjustments to your bike. Your special need is a challenge for us!

Accessories & options

The Copilot 3 can be expanded with various accessories and options. Think of a another seat, handlebars or pedals. View the options on our website or inquire with your supplier.

Colors

The tandem Copilot 3-24 (supervisor with child in front,

from 4 years of age) is available in white (RAL9016) and pearl dark grey (RAL9023).

The tandem Copilot 3-26 (adult supervisor in front) is available in white (RAL9016) and pearl dark grey (RAL9023).

You want to design your Copilot 3 even more personal? We can make the Copilot 3 in any desired RAL color for you!



The most important specifications

Technical specifications	Copilot 3-24		Copilot 3-26	
	front	rear	front	rear
Length in mm	2040-2080		2150	
Width in mm	890		890	
Weight in kg	30		31	
Crank length in mm	150	170	150	170
Length inner leg in mm	500-780	710-900	690-970	850-990
Seat height to ground in mm	640-780	730-940	860-1000	940-1080
Maximum load in kg	180		180	
Tire size front in ETRTO	47-406	47-507	47-406	47-559
Gears	8-speed hub with coaster brake (Shimano Nexus)			

About Huka

Developed in Oldenzaal, we have been producing bikes that make a difference in the world of riders for more than 40 years. Your mobility, independence and freedom is our passion. With a committed team, we provide a wide collection of innovative bicycles and other means of transport with Dutch design. You want to learn more about the Copilot 3 and our assortment? Give us a call or visit our website!