

The Huka logo is a yellow diamond shape with the word 'huka' written in a dark grey, lowercase, sans-serif font.A yellow circular badge with the text 'TAKE A FREE TESTRIDE!' in a bold, black, uppercase, sans-serif font.

TANDEM COPILOT

## Suitable for children & adults

Cycle in front, for the best cycling experience!

Are you looking for a tandem where the passenger has the best cycle experience and the attendant has the most control? Experience the Copilot tandem! It is the ideal solution if you still want to cycle, but when this is no longer possible without help. The Copilot gives both cyclists the freedom and control you want. The tandem can also be adjusted to specific preferences and wishes. With the Copilot tandem you can cycle together without any problems!

We look at individual cyclists and their needs and want to match them as well as possible. This way, we continuously innovate and optimize our bikes, so that your bike is equipped with the best materials and optimum comfort. For example, the attendant is positioned higher than the passenger to keep a good overview of the traffic. In addition, you can easily get on and off with the Copilot. The lighting ensures that you can see well and are being seen in the dark.

The tandem Copilot is equipped with a low step-in entry, easy-to-use hand brakes, a large closed modern chain guard and an AXA cable lock.

- ✔ Available with electrical support.
- ✔ Both the co-steering function and the co-peddalling function are adjustable for the passenger.
- ✔ Available in two models: supervisor adults & supervisor child.
- ✔ Supervisor cycles on the back for optimum view of the passenger.
- ✔ 10-year warranty on the frame and the front.\*

\* Register your Huka bike once via [www.huka.nl/en/warranty](http://www.huka.nl/en/warranty) and get a 10 year warranty on the frame and front fork. Standard you get a 5 year warranty on these parts and 2 years on other components (with the exception of wear sensitive parts).

### Electrical support

The Copilot is available with an auxiliary motor (PAS vario). With the help of this drive it is easier to pedal and therefore less tiring. This silent pedal support is standard supplied with a drive help, three travel speed settings and a clear display.

### Support & special adjustments

In addition to the extra accessories that we can offer, such as hip support, trunk support, foot fixations etc.

We also provide customized solutions. That way we can make specific adjustments to your bike. Your special need is a challenge for us!

### Accessories & options

The Copilot can be expanded with various accessories and options. Think of another seat, handlebars or pedals. View the options on our website or inquire with your supplier.

### Colours

The tandem Copilot 24 (attendant with child in front, from 4 years of age) is available in white (RAL9016) and pearl dark grey (RAL9023).

The tandem Copilot 26 (adult supervisor in front) is available in white (RAL9016), pearl dark grey (RAL9023).

Any other RAL colours are available on request.



### The most important specifications

Technical specifications	Copilot 24		Copilot 26	
	front	rear	front	rear
Length in mm	2470		2580	
Width in mm	550	630	650	650
Weight in kg	29		31	
Crank length in mm	150	170	170	
Length inner leg in mm	410-790	710-940	560-960	770-930
Seat height to ground in mm	640-770	810-990	810-1020	900-1060
Maximum load in kg	180		180	
Tire size front in ETRTO	47-507	47-507	47-559	47-559
Gears	8-speed hub with coaster brake (Shimano Nexus)			

### About Huka

Developed in Oldenzaal, we have been producing bikes that make a difference in the world of riders for more than 40 years. Your mobility, independence and freedom is our passion. With a committed team, we provide a wide collection of innovative bicycles and other means of transport with Dutch design. You want to learn more about the Copilot and our assortment? Give us a call or visit our website!

