



**huka**

TAKE A  
FREE  
TESTRIDE!

**WHEELCHAIR BIKE DUET**

## Bicycle & wheelchair in one

### Going out together in the open air

The wheelchair bicycle Duet is a unique combination of a bicycle and wheelchair that we have developed for both adults and children. One cycles and the other sits in the wheelchair at the front. Thanks to the design, pedalling and steering is light and smooth. Arrived at your destination? Simply disconnect the chair and it will become a separate manual wheelchair. Do you have specific requirements for the wheelchair bicycle? We adjust the wheelchair bicycle entirely to your circumstances. It is a smart solution in all kinds of daily situations!

We look at individual cyclists and their needs and want to match them as well as possible. This way, we continuously innovate and optimize our bikes, so that your bike is equipped with the best materials and optimum comfort. For example, the wheelchair has stable driving characteristics, due to the parking brake, and it is easy to get on and off.

The Duet wheelchair bicycle is standard equipped with, among other things, 3 gears, double drum brakes, coaster brake, parking brake, hip belt for the wheelchair and seat cushions.

- ✔ Available with electrical support the wheelchair part can be disconnected and can be used as a (push) wheelchair.
- ✔ Thanks to various adjustment options, the bicycle can be adjusted to the posture of the passenger.
- ✔ 10-year warranty on the frame.\*

\* Register your Huka bike once via [www.huka.nl/en/warranty](http://www.huka.nl/en/warranty) and get a 10 year warranty on the frame. Standard you get a 5 year warranty on this part and 2 years on other components (with the exception of wear sensitive parts).

## Electrical support

The Duet is available with an auxiliary motor (PASvario). With the help of this drive it is easier to pedal and therefore less tiring. This silent pedal support is standard supplied with a drive help, three travel speed settings and a clear display.

## Support & special adjustments

In addition to the extra support that we can offer, such as side cushions, with which the seat depth and width can be adjusted,

we also provide customized solutions. The chair is modified to the user (adults or child). Do you have a special custom-made seat? Then the seat shell can be replaced by the orthosis. Your special need is a challenge for us!

## Accessories & options

The Duet can be expanded with various accessories and options. Think of a headrest, spokes cover, four-point belt or seat pants. That is no problem for us!


View the options on our website or inquire with your dealer.

## Colors

The Duet is available in the following colours:

- Silver grey
- Metallic red
- Metallic blue

The front seat is supplied in yellow with blue insert cushions. Any other RAL colours are available on request.



## The 10 most important specifications

Technical specifications	bicycle part	wheelchair part
Crank length in mm	170	n/a
Length inner leg in mm	780-930	100-550
Tire size in inch/ ERTRO	26 x 1,75 (47-559)	26 x 1,75 (47-559)
Maximum load in kg	100	100
Total length x width in mm	2350 - 650	
Maximum load total in kg	180	
Seat width x seat depth in mm	n/a	450 x 400
Length inner leg in mm	n/a	100 - 550
Back height in mm	n/a	550
Maximum load in kg	100	100

## About Huka

Developed in Oldenzaal, we have been producing bikes that make a difference in the world of riders for more than 40 years. Your mobility, independence and freedom is our passion. With a committed team, we provide a wide collection of innovative bicycles and other means of transport with Dutch design. You want to learn more about the Duet and our assortment? Give us a call or visit our website!