

Pricelist/ order form Duet



Name _____

Organisation _____

Street/ Nr. _____

Postal code/ City _____

Telephone _____

Base Duet

The Duet basically comes with a number of standard parts. These parts are marked with a symbol in front of this specific part. The Duet will be provided with these parts if no other option from this category has been selected. Also a number of parts come with a standard addition. These parts are marked with a symbol behind. It will have this standard addition if no alternative has been selected. For the complete description of the standard parts we refer to the specifications. Below you will also find all possible options and accessories.

	Art.nr.	Base	Price €
S	450000	Wheelchair bicycle Duet	4.643
	450010	- Frame color metallic blue	-
	450011	- Frame color metallic red	-
	450012	- Frame color silver grey	-

	Art.nr.	Drives (not applicable with electric drive)	Price €
S	490081	3 speed hub with pedal brake	-
	490080	8-speed hub Shimano Nexus with pedal brake	287
	490082	8-speed hub Shimano Nexus with pedal brake automatic	512
	490083	Rohloff 14 speed gear hub freewheel	1.875

	Art.nr.	Drives (not applicable with electrical support)	Price €
	408175	Mountaindrive system in bottom bracket	753

	Art.nr.	(Electrical) support	Price €
S	469800	Without electrical support	-
	469801	Pedal support PAS-vario direct drive, 36 Volt (incl. 14,25 Ah Li-Ion accu), 7 speed	2.693
	400063	S - Without pedal brake for Duet electro	-
	400063	Pedal brake for bicycle with electrical support	528

	Art.nr.	Handlebar attachments	Price €
S	408006	Handlebar bracket, standard	-
	408012	Handlebar bracket, 20 cm extended	84

	Art.nr.	Seats	Price €
S	401637	Velo Plush seat	-
	402860	Comfort gel seat, unisex	62
	426837	Female seat, leather version	158
	426836	Male seat, leather version	158
	418265	Wide seat	62
	426830	Skirt seat	44

Make the Wheelchair bicycle Duet even more personal? Special adjustments are possible on request!

The mentioned prices are exclusive VAT, delivery ex works. The prices are valid from July 1, 2021 and replace all previous price lists. Subject to technical changes and printing errors.

Pricelist/ order form Duet

	Art.nr.	Seat posts	Price €
S	432059	Suspension seat post, length 350 mm	-
	409033	Aluminium seat post, length 350 mm	-

	Art.nr.	Seating system wheelchair	Price €
S	491000	Seat shell with belt and foot support, color yellow	-
	491063	Without seat shell belt, foot support and frame foot support	-270
	491064	Mounting plate incl. adjustable adapter	522

	Art.nr.	Wheelchair seat adjustments	Price €
	491001	Adjustable headrest with cushions, color yellow	218
	491011	Four-point safety belt	79
	491012	Breeches	105
	491013	Foot strap, pair	41
	491014	Calf band	41
	491009	Footrest raise (for children or very long legs)	170
	491041	Back cushion for reducing seat depth by 5 cm	142
	491040	Back cushion for reducing seat depth by 10 cm	142
	491042	Side supports for reducing seat width by 8 cm	175
	491043	Side supports for reducing seat width by 12 cm	175
	491044	Side supports for reducing seat width by 20 cm	175
	491060	Backrest shell reinforced and fixed angle	84
	491006	Raised armrest, upholstered and removable	169

	Art.nr.	Other adjustments	Price €
	491003	Spoke covers wheelchair part (unbreakable and transparent), per pair	248
	491094	Extra parking brakes placed on wheelchair	292
	491007	Hand protection	185
	491023	Quick release footboard	159
	491024	Quick release handlebar	79
	491025	Quick release headrest	40
	426961	Rearview mirror	35
	491046	Anti-tip wheel, connectable	217
	400055	Cycle computer (display speed, day counter and mileage)	85
	430001	Wheelchair poncho	138

Date of order _____

Signature/ stamp

Total €	
VAT	
Total € incl. VAT	

Make the Wheelchair bicycle Duet even more personal? Special adjustments are possible on request!

The mentioned prices are exclusive VAT, delivery ex works. The prices are valid from July 1, 2021 and replace all previous price lists. Subject to technical changes and printing errors.