




## Technical data Pendel

	Standard version
Material chassis	tubular steel frame
Housing	polyester material
Colors <sup>A</sup>	available in 6 colors (see price list)
Front wheel suspension	with suspension
Rear wheels suspension	gas suspension/ air suspension
Handlebar <sup>A</sup>	steel, height adjustable
Front brakes <sup>A</sup>	mechanical drum brake
Rear brakes <sup>A</sup>	hydraulic drum brake
Electronics	Sevcon, programmable
Batteries	2 x 12 V, 105 Ah
Charger	24 V, 12 Ah
Motor	2200 Watt
Options	see price list for different options
Guarantee	2 years, with the exception of parts that will wear

	Standard version
Total length in mm	1760
Length inside in mm	1140
Width total in mm (standard / wide version)	1000/ 1080
Maximum wheelchair width in mm (std / wide version)	720/ 800
Ground clearance in mm	110
Reversal space in mm (standard / wide version)	2860/ 2900
Turning circle in mm	3900
Maximum rise in degrees	17
Maximum load in kg	125
Total weight Pendel incl. batteries in kg	160
Tire size front	3.00 -10
Tire size rear	3.00 -12
Maximum speed as a medical device in km/h	15
Maximum speed as a vehicle in km/h	25
Range in km	40-60

We continuously innovate and optimize our bicycles and other means of transport, so that your Huka is equipped with the best materials and optimum comfort. As a result, the delivered product may differ from the images shown. The accessories shown are not included. Printing errors reserved.