




Technical data City, size M

	Standard version
Material frame	anti-corrosion treated steel
Coating	outside quality epoxy powder coating
Color A	light grey (RAL7047) or azure blue (RAL5009)
Wheels	aluminum rim and galvanized hub with stainless steel spokes
Tires	Schwalbe Marathon Plus (anti-puncture)
Lighting	LED battery lighting front and rear
Bell A	sporty and easy-to-use bell, type Decibel
Mudguards	impact and form-retaining plastic
Chain protection	Hebie Chainglider
Handlebar and stem A	slightly bent handlebar, ergonomic stem
Handles	Ergo Grip
Seat A	Velo Plus, ergonomic comfort seat
Seatpost A	aluminium, 350 mm, height adjustable
Brakes front	Shimano v-brake
Brakes rear	hydraulic disc brakes incl. parking brake
Drive A	freewheel with differential in rear axle
Key lock	AXA Solid security lock
Options	see price list for different options

	Standard version
Total length in mm (H)	1675
Width total in mm	785
Weight in kg	25
Crank length in mm (A)	170
Length inner leg in mm (B)* **	660-880
Length inner leg with lowered seatpost in mm (B)**	570-660
Seat height to ground in mm (C)	min. 720
Height of luggage carrier (G)	610
Entry space in mm (D)	570
Entry height in mm (E)	190
Height adjustment handlebar in mm (F)	110
Handlebar width in mm	640
Maximum load in kg	130
Tire size front in ETRTO	47-507
Tire size front in inches	24 x 1,75
Tire size rear in ETRTO	47-406
Tire size rear in inches	20 x 1,75

* When mounting the backrest, the inner leg length increased by 50 mm in the lowest position
 **Adjustable in combination with crank shorteners



Warranty

Huka has a standard 5-year warranty on the frame and front fork. If you register your bike via www.huka.nl/en/warranty, you are eligible to receive a 10-year warranty on these parts. Additionally, you have a 2-year warranty on other parts with the exception of wear-sensitive parts. You can find more information about the City warranty on www.huka.nl/en/warranty.

Possibilities for alternative parts

The City is equipped with a number of standard parts. Optional these parts can be replaced by alternative parts. The relevant parts are marked by a symbol behind this specific part. The alternative parts can be found on the price list.

The delivered product may differ from the images shown. The accessories shown are not included. Printing errors reserved.