

hukaTAKE A  
FREE  
TESTRIDE!

CUSTOM CHILDREN AT-SERIES

Happy & safe for  
children

## Specially developed tricycle for children

Are you looking for a children's bike that offers more stability and support? With the adapted AT bike your child can enjoy cycling without difficulties. Playing outside, cycling to school and playing with friends? Easy! This children's bicycle can also be adjusted entirely to the wishes and circumstances of your child. For example, the lightweight bicycle can be supplied in all hip, cool and cheerful RAL colours and there are various accessories with which you can make the bicycle even more personal and attractive. Your child can cycle carefree!

We look at individual cyclists and their needs and want to match them as well as possible. This way, we continuously innovate and optimize our bikes, so that your bike is equipped with the best materials, latest technology, and optimum comfort. We have two types that have been specially developed for two age groups, the AT-A (up to and including about 5 years) and the AT-B bike (from about 5 years up to and including 9 years).

The AT bike is standard equipped with, among other things, a closed modern chain guard, lightly operated Rollerbrake hand brakes and a low and wide step-on plate. The AT-B is also equipped with lighting and a safety lock (AXA solid).

- ✓ Available in two sizes: AT-A and AT-B.
- ✓ Safe to get on and off and easy to maneuver, due to its three-wheel construction and low center of gravity.
- ✓ The light running bike with stable handling. Individually adaptable with various supports and accessories.
- ✓ Cheerful and child-friendly styled.
- ✓ 10 year warranty on the frame and front fork.\*

\*After registration at [www.huka.nl/en/warranty](http://www.huka.nl/en/warranty).



## Electrical support

The ATB is available with an auxiliary motor (PASvario). With the help of this drive it is easier to pedal and therefore less tiring. This silent pedal support is supplied with five travel speed settings and a well-arranged display.

## Make your AT bike personal

In addition to the extra support that we can offer, such as a push bar, hip support, trunk support, a back support, foot

fixations, etc. We also provide customized solutions. This allows us to make personal adjustments to the children's bicycle. Special needs are a challenge for us!

## Accessories & options

The AT can be equipped with various accessories and options. Think of a different seat, handlebars or pedals. View the options on our website or contact your supplier.

## Colours

The AT-A tricycle is standard supplied in light blue (RAL5024). The AT-B tricycle standard comes in the following colours:

- Bright blue (RAL5015)
- Pink (RAL3014)

Want to make the AT even more personal with, for example, your child's favourite colours? That is possible! We can make the AT children's bike in any hip, cool or cheerful RAL colour.



## The 10 most important specifications

Technical specifications	AT-A	AT-B
Length in mm	1140	1360
Width in mm	540	625
Weight in kg	17	24
Crank length in mm	90	110
Length inner leg in mm	380-490	470-650
Entry frame length in mm	300	340
Entry height frame in mm	145	260
Maximum load in kg	40	60
Tire size front in inch / ETRTO	16 x 1,75/ 47-305	20x 1,75/ 35-406
Tire size rear in inch / ETRTO	12 x 1,95/ 50-203	16 x 1,75/ 47-305

## About Huka

Developed in Oldenzaal, we have been producing bikes that make a difference in the world of riders for almost 50 years. Your mobility, independence and freedom are our passion. With a committed team, we provide a wide collection of innovative bicycles Dutch design. You want to learn more about the AT bike and our assortment? Give us a call or visit our website!