

The Huka logo is displayed in a yellow, tilted rectangular box in the top left corner. The word "huka" is written in a bold, lowercase, sans-serif font.A yellow circular badge with the text "TAKE A FREE TESTRIDE!" in a bold, uppercase, sans-serif font, located in the top right corner of the image.

**TANDEM COPILOT 3**

# Stable getting on and off with three wheels

## Extra safe and stable with the passenger in front

Are you looking for a tandem that is extra stable and safe for, for example, a restless mobile passenger or a passenger with a reduced sense of balance? Experience the tandem Copilot 3. With the two wheels at the front you have the right balance and support which is needed for this passenger. Suitable for adults and children! Thanks to the third wheel, both riders can remain seated even when the tandem is standing still. With this tandem, only the supervisor can steer, and the passenger has a handlebar for guidance. The Copilot 3 can also be personalized to specific preferences and wishes. It is a tandem to enjoy cycling together without worrying about safety!

We look at individual cyclists and their needs and want to match them as well as possible. This way, we continuously innovate and optimize our bikes, so that your bike is equipped with the best materials, latest technology and optimum comfort. For example, the supervisor is positioned higher than the passenger to keep a good overview of the traffic. In addition, you can easily get on and off with the Copilot 3. The supervisor on the back can optionally decide whether the passenger has to pedal along or not.

The tandem Copilot 3 is standard equipped with a low step-up entry, easy-to-use hand brakes, a large closed modern chain guard and an AXA cable lock.

- ✓ Extremely stable due to the two wheels in front. Also for mobile passengers.
- ✓ Available in two models: supervisor - adult & supervisor- child.
- ✓ Standard equipped with a removable and tiltable front frame, making it easier to get on and store the bike more compactly
- ✓ 10-year warranty on the frame and the front fork.\*

\*After registration at [www.huka.nl/en/warranty](http://www.huka.nl/en/warranty).



### Electrical support

The Copilot 3 is available with an auxiliary motor (PASvario). With the help of this drive it is easier to pedal and therefore less tiring. This silent pedal support is standard supplied with a start help, five travel speed settings and a clear display.

### Make you Copilot 3 personal

In addition to the extra accessories that we can offer, such as hip support, trunk support, foot fixations, etc. We also provide customized

solutions. That way we can make personal adjustments to your bike. Your special need is a challenge for us! The Copilot 3 can be expanded with various accessories and options. Think of a another seat, handlebars or pedals. View the options on our website or contact your supplier.

### Frame at front

The frame at the front is detachable and can be replaced by a smaller or a larger part. In this way, the Copilot 3, can be used flexibly and can also be

used as a grow along bike.

### Colors

The tandem Copilot 3-24 (co-driver in front, child - starting at the age of 4 years) and the tandem Copilot 3-26 (adult supervisor in front) are is available in white (RAL9016) and pearl dark grey (RAL9023).

You want to design your Copilot 3 even more personal? We can make the Copilot 3 in any desired RAL color for you!



### The most important specifications

Technical specifications	Copilot 3-24		Copilot 3-26	
	front	rear	front	rear
Length in mm	2040-2080		2150	
Width in mm	890		890	
Weight in kg	30		31	
Crank length in mm	150	170	150	170
Length inner leg in mm	500-780	710-900	690-970	850-990
Seat height to ground in mm	640-780	730-940	860-1000	940-1080
Maximum load in kg	180		180	
Tire size front in ETRTO	47-406	47-507	47-406	47-559
Gears	8-speed hub with coaster brake (Shimano Nexus)			

### About Huka

Developed in Oldenzaal, we have been producing bikes that make a difference in the world of riders for almost 50 years. Your mobility, independence and freedom are our passion. With a committed team, we provide a wide collection of innovative bicycles with Dutch design. You want to learn more about the Copilot 3 and our assortment? Give us a call or visit our website!