

STABLE GETTING ON AND OFF WITH THREE WHEELS

# Tandem Copilot 3

## Extra safe and stable with the passenger in front

*Are you looking for a tandem that is extra stable and safe for, for example, a restless mobile passenger or a passenger with a reduced sense of balance? Experience the tandem Copilot 3. With the two wheels at the front you have the right balance and support which is needed for this passenger. Suitable for adults and children! Thanks to the third wheel, both riders can remain seated even when the tandem is standing still. With this tandem, only the supervisor can steer, and the passenger has a handlebar as a handhold. The Copilot 3 can also be personalized to specific preferences and wishes. It is a tandem to enjoy cycling together without worrying about safety!*

We look at individual cyclists and their needs and want to match them as well as possible. This way, we continuously innovate and optimize our bikes, so that your bike is equipped with the best materials, latest technology and optimum comfort. For example, the supervisor is positioned higher than the passenger to keep a good overview of the traffic. In addition, you can easily get on and off with the Copilot 3. The supervisor on the back can optionally decide whether the passenger has to pedal along or not. The tandem Copilot 3 is standard equipped with a low step-up entry, hydraulic brakes at the front, a large closed modern chain guard and an AXA cable lock.

- ✓ Extremely stable due to the two front wheels, allowing for extra safe cycling with a mobile passenger or a passenger with a reduced sense of balance.
- ✓ Available with and without electrical assistance (Bafang) in two models: supervisor - child (24") & supervisor - adult (26")
- ✓ Standard equipped with a removable and tiltable front frame, making it easier to get on and store the bike more compactly
- ✓ Equipped with rear hydraulic disc brakes, which require less hand power.
- ✓ 10-year warranty on the frame and the front. After registration at [www.hukabikes.com/warranty](http://www.hukabikes.com/warranty).

### Electrical support

The Copilot 3 is available with an auxiliary motor (PASvario). With the help of this drive it is easier to pedal and therefore less tiring. This silent pedal support is standard supplied with a start help, five travel speed settings and a clear display.

### Make you Copilot 3 personal

In addition to the extra accessories that we can offer, such as hip support, trunk support, foot fixations, etc. We also provide customized solutions. That way we can make personal adjustments to your bike. Your special need is a challenge for us!

The Copilot 3 can be expanded with various accessories and options. Think of a another seat, handlebars or pedals. View the options on our website or contact your supplier.

### Frame at front

The frame at the front is detachable and can be replaced by a smaller or a larger part. In this way, the Copilot 3, can be used flexibly and can also be used as a grow along bike.

### Colors

The tandem Copilot 3-24 (co-driver in front, child - starting at the age of 4 years) and the tandem Copilot 3-26 (adult supervisor in front) are available in white (RAL9016) and pearl dark grey (RAL9023). You want to design your Copilot 3 even more personal? We can make the Copilot 3 in any desired RAL color for you!



## The most important specifications

Technical specifications	Copilot 3-24		Copilot 3-26	
	front	rear	front	rear
Material frame & coating	steel, epoxy coated			
Electric support	option, Bafang rear wheel motor 65 Nm			
Tyre size front in inch (ETRTO)	20 x 1,75 (47-406)		20 x 1,75 (47-406)	
Tyre size rear in inch (ETRTO)	24 x 1,75 (47-507)		26 x 1,75 (47-559)	
Inner leg length in mm	500-780	710-900	690-970	850-990
Step up height in mm	400		400	
Saddle height to ground in mm	620-870	770-940	860-1000	940-1080
Maximum load in kg	180		200	
Length x width in mm	2040/2080 x 890		2150 x 890	

### About Huka

Developed in Oldenzaal, we have been producing bikes that make a difference in the world of riders for almost 50 years. Your mobility, independence and freedom are our passion. With a committed team, we provide a wide collection of innovative bicycles with Dutch design. You want to learn more about the Copilot 3 and our assortment? Give us a call or visit our website!